

Product datasheet

Recombinant Human Mago nashi homolog 2 protein ab156754

[1 Image](#)

Description

Product name	Recombinant Human Mago nashi homolog 2 protein
Purity	> 90 % SDS-PAGE. ab156754 was purified using conventional chromatography techniques.
Expression system	Escherichia coli
Accession	Q96A72
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	MGSSHHHHHH SSGLVPRGSH MGSMASVDF YLRYYVGHKG KFGHEFLEFE FRPDGKLRYA NNSNYKNDVM IRKEAYVHKS VMEELKRIID DSEITKEDDA LWPPPDRVGR QELEIVIGDE HISFTTSKIG SLIDVNQSKD PEGLRVFYLL VQDLKCLVFS LIGLHFKIKP I
Predicted molecular weight	20 kDa
Amino acids	1 to 148
Tags	His tag N-Terminus

Specifications

Our [Abpromise guarantee](#) covers the use of **ab156754** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	Mass Spectrometry SDS-PAGE
Mass spectrometry	MALDI-TOF
Form	Liquid

Preparation and Storage

Stability and Storage

Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.


pH: 8.00

Constituents: 0.02% DTT, 0.32% Tris HCl, 30% Glycerol (glycerin, glycerine), 0.88% Sodium chloride

General Info

Function	Involved in mRNA splicing and in the nonsense-mediated decay (NMD) pathway.
Tissue specificity	Ubiquitous.
Sequence similarities	Belongs to the mago nashi family.
Cellular localization	Nucleus.

Images



15% SDS-PAGE analysis of ab156754 (3µg).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors