

Recombinant Human MS4A4A protein ab162515

1 Image

Description	
Product name	Recombinant Human MS4A4A protein
Expression system	Wheat germ
Accession	<u>Q96JQ5</u>
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	MTTMQGMEQAMPGAGPGVPQLGNMAVIHSHLWKGLQEK FLKGEPKVLGVV QILTALMSLSMGITMMCMASNTYGSNPISVYIGYTWGSVMF IISGSLSI AAGIRTTKGLVRGSLGMNITSSVLAASGILINTFSLAFYSFHH PYCNYYG NSNNCHGTMSILMGLDGMVLLLSVLEFCIAVSLSAFGCKV LCCTPGGVVL ILPSHSHMAETASPTPLNEV
Amino acids	1 to 220
Tags	GST tag N-Terminus

Specifications	
Our <u>Abpromise guarantee</u> covers the use of ab162515 in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
Applications	Western blot ELISA
Form	Liquid
Additional notes	

Preparation and Storage	
Stability and Storage	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.

pH: 8.00

Constituents: 0.31% Glutathione, 0.79% Tris HCl

General Info


Relevance

MS4A4A is a member of the membrane spanning 4A gene family. This nascent protein family are characterized by common structural features and similar intron/exon splice boundaries and display unique expression patterns among hematopoietic cells and non lymphoid tissues. MS4A4A may be involved in signal transduction as a component of a multimeric receptor complex. There are two named isoforms.

Cellular localization

Cell Membrane; Multi pass membrane protein.

Images



SDS-PAGE - Recombinant Human MS4A4A protein
(ab162515)

ab162515 on a 12.5% SDS-PAGE stained with Coomassie Blue.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors