

Product datasheet

Recombinant Human NANP protein ab99392

1 Image

Description

Product name	Recombinant Human NANP protein
Purity	> 90 % SDS-PAGE. ab99392 is purified using conventional chromatography techniques.
Expression system	Escherichia coli
Accession	<u>Q8TBE9</u>
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	MRGSHHHHHHGMASMTGGQQMGRDLYDDDDKDRWG S MGLSRVRAVFFDLNLTIDTAGASRRGMLEVIKLLQSKYHY KEEAEIICD KVQVKLSKECFHPYNTCITDLRTSHWEEAIQETKGAANR KLAEECYFLW KSTRLQHMTLAEDVKAMLTELKRVRLLLL TNGDRQTQRE KIEACACQSY FDAVVVGGEQREEKPAPSIYCCNLLGVQP GDCVMVG DTLETDIQGGLN AGLKATVWINKNGVPLKSSPVPHYMVSSVLELPALLQSID CKVSMST
Predicted molecular weight	32 kDa including tags
Amino acids	1 to 248
Tags	His tag N-Terminus

Specifications

Our **Abpromise guarantee** covers the use of **ab99392** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	SDS-PAGE Mass Spectrometry
Mass spectrometry	MALDI-TOF

Form Liquid

Preparation and Storage


Stability and Storage Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
pH: 8.00
Constituents: 0.0308% DTT, 0.316% Tris HCl, 10% Glycerol (glycerin, glycerine), 0.58% Sodium chloride

General Info

Pathway Amino-sugar metabolism; N-acetylneuraminate biosynthesis.

Sequence similarities Belongs to the HAD-like hydrolase superfamily. NANP family.

Images



15% SDS-PAGE analysis of 3µg ab99392.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors