

## Product datasheet

# Recombinant Human NAPSIN A protein (His tag) ab219242

[1 Image](#)

### Description

<b>Product name</b>	Recombinant Human NAPSIN A protein (His tag)
<b>Purity</b>	> 90 % SDS-PAGE. ab219242 was purified by conventional chromatography.
<b>Expression system</b>	Escherichia coli
<b>Accession</b>	<u><b>O96009</b></u>
<b>Protein length</b>	Full length protein
<b>Animal free</b>	No
<b>Nature</b>	Recombinant
<b>Species</b>	Human
<b>Predicted molecular weight</b>	41 kDa including tags
<b>Amino acids</b>	64 to 420
<b>Tags</b>	His tag N-Terminus
<b>Additional sequence information</b>	This is the mature protein without the signal sequence or propeptide sequence (NP_004842).

### Specifications

Our **Abpromise guarantee** covers the use of **ab219242** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

<b>Applications</b>	SDS-PAGE Mass Spectrometry
<b>Mass spectrometry</b>	MALDI-TOF
<b>Form</b>	Liquid

### Preparation and Storage

<b>Stability and Storage</b>	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle. pH: 8.00 Constituents: 0.32% Tris HCl, 10% Glycerol (glycerin, glycerine)
------------------------------	--

## General Info

---

<b>Function</b>	May be involved in processing of pneumocyte surfactant precursors.
<b>Tissue specificity</b>	Expressed predominantly in adult lung (type II pneumocytes) and kidney and in fetal lung. Low levels in adult spleen and very low levels in peripheral blood leukocytes.
<b>Sequence similarities</b>	Belongs to the peptidase A1 family.
<b>Cellular localization</b>	Secreted.

---

## Images

---



15% SDS-PAGE analysis of ab219242 (3µg).

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors