

Recombinant Human NAT13 protein ab123145

1 Image

Description

Product name	Recombinant Human NAT13 protein		
Purity	> 95 % SDS-PAGE. ab123145 was purified using conventional chromatography.		
Expression system	Escherichia coli		
Accession	<u>Q9GZZ1</u>		
Protein length	Full length protein		
Animal free	No		
Nature	Recombinant		
Species	Human		
Sequence	MGSSHHHHHHSSGLVPRGSHMGSHMKGSRIELGDVTPH NIKQLKRLNQVI FPVSYNDKFYKDVLEVGE LAKLAYFN DIAVGAVCCRVDHS QNQKRLYIMT LGCLAPYRRLGIGTKMLNHVLNICEKDGTFDNMLHVQISNE SAIDFYRK FGFEIIE TKKNYYKRIEPADAHVLQKNLKVPSGQNADVQKT DN		
Predicted molecular weight	22 kDa including tags		
Amino acids	1 to 169		
Tags	His tag N-Terminus		

Specifications

Our **Abpromise guarantee** covers the use of **ab123145** in the following tested applications.
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	SDS-PAGE Mass Spectrometry
Mass spectrometry	MALDI-TOF
Form	Liquid

Preparation and Storage

Stability and Storage

Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 8.00

Constituents: 0.02% DTT, 0.32% Tris HCl, 10% Glycerol (glycerin, glycerine), 0.58% Sodium chloride

General Info

Function

Probable catalytic component of the NAA11-NAA15 complex which displays alpha (N-terminal) acetyltransferase activity.


Sequence similarities

Belongs to the acetyltransferase family. GNAT subfamily.
Contains 1 N-acetyltransferase domain.

Cellular localization

Cytoplasm.

Images



15% SDS PAGE, 3 µg of ab123145 loaded.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors