

Recombinant Human NIP7 protein ab95504

1 Image

Description

Product name	Recombinant Human NIP7 protein
Purity	> 90 % SDS-PAGE. ab95504 is purified using conventional chromatography techniques.
Expression system	Escherichia coli
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	MRPLTEEETR VMFEKIAKYI GENLQLLVDR PDGTYCFRLH NDRVYYVSEK IMKLAANISG DKLVSLGTCFGKFTKTHKFR LHVTALDYL PYAKYKWWIK PGAEQSFLYG NHVLKSGLGR ITENTSQYQG VVVYSMADIPLGFGVAAKST QDCRKVDPM IIVFHQADIG EYVRHEETLT LEHHHHHH

Specifications

Our **Abpromise guarantee** covers the use of **ab95504** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	SDS-PAGE Mass Spectrometry
Form	Liquid


Preparation and Storage

Stability and Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. pH: 8.00 Constituents: 0.316% Tris HCl, 20% Glycerol (glycerin, glycerine), 0.58% Sodium chloride
------------------------------	---

General Info

Function	Required for proper 34S pre-rRNA processing and 60S ribosome subunit assembly.
Tissue specificity	Expressed in hematopoietic stem/progenitor cells.
Sequence similarities	Belongs to the NIP7 family. Contains 1 PUA domain.
Cellular localization	Nucleus > nucleolus.

Images



15% SDS-PAGE showing ab95504 at approximately 21.5kDa (3µg).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors