

Product datasheet

Recombinant Human NuMA protein ab114324

1 Image

Description

Product name	Recombinant Human NuMA protein
Expression system	Wheat germ
Accession	<u>Q14980</u>
Protein length	Protein fragment
Animal free	No
Nature	Recombinant
Species	Human
Sequence	SPASPMGDILQTPQFQMRRLKKQLADERSNRDELELELA ENRKLLTEKDA QIAMMQQRIDRLALLNEKQAASPLEPKEELELRDKNESLT MRLHETLKQC QDLKTEK
Predicted molecular weight	37 kDa
Amino acids	200 to 306

Specifications

Our **Abpromise guarantee** covers the use of **ab114324** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	SDS-PAGE
	Western blot
	ELISA
Form	Liquid

Preparation and Storage

Stability and Storage	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.
	pH: 8.00
	Constituents: 0.3% Glutathione, 0.79% Tris HCl


General Info

Function

May be a structural component of the nucleus.

Cellular localization

Nucleus. Chromosome. Dissociates from condensing chromosomes during early prophase, before the complete disintegration of the nuclear lamina. As mitosis progresses it reassociates with telophase chromosomes very early during nuclear reformation, before substantial accumulation of lamins on chromosomal surfaces is evident.

Images

SDS-PAGE analysis of ab114324 on a 12.5% gel stained with Coomassie Blue.

SDS-PAGE - Recombinant Human NuMA protein
(ab114324)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors