

Recombinant Human PCQAP / MED15 protein  
ab162600

1 Image

Description	
Product name	Recombinant Human PCQAP / MED15 protein
Expression system	Wheat germ
Accession	<u>Q96RN5</u>
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	MDVSGQETDWRSTAFRQKLVSQIEDAMRKAGVAHSKSS KDMESHVFLKAK TRDEYLSLVARLIHFRDIHNKKSQASVSAQLQLQQVALQQ QQQQQQQFQQ QQQAALQQQQQQQQQFQAQQSAMQQQFQAVVQQQQQ LQQQQQQQHLLKLH HQNQQQIQQQQQQLQRIQLQLQQQQQQQQQQQQQQQ ALQAQPPIQQPP MQQPQPPPSQALPQQLQQMHHHTQHHQPPPQPQPPVA QNQPSQLPPQSQT QPLVSQAQALPGQMLYTQPPLKFVRAPMVVQPPVQPQ VQQQQTAVQTAQ AAQMVAPGVQVSQSSLPMLSSPSPGQQVQTPQSMPPP PQPSPQPGQPSSQ PNSNVSSGPAPSPSSFLPSPSPQPSQSPVTARTPQNFSV PSPGPLNTPVN PSSVMSPAGSSQAEEQQYLDKLKQLSKYIEPLRRMINKID KNEDRKKDLS KMKSLLDILTDPSKRCPLKTLQKCEIALEKLKNDMAVPTTP PPVPPTKQ QYLCQPLLDVLANIRSPVFNHSLYRTFVPAMTAIHGPPITA PVVCTRKR RLEDDERQSIPSVLQGEVARLDPKFLVNLDPSHCSNNGT VHLICKLDDKD LPSVPPLELSVPADYPAQSPLWIDRQWQYDANPFLQSVH RCMTSRLLQLP DKHSVTALLNTWAQSVHQACLSAA

<b>Amino acids</b>	1 to 674
<b>Tags</b>	GST tag N-Terminus

## Specifications

Our **Abpromise guarantee** covers the use of **ab162600** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

<b>Applications</b>	Western blot
	ELISA
<b>Form</b>	Liquid

## Additional notes

## Preparation and Storage


<b>Stability and Storage</b>	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.
	pH: 8.00
	Constituents: 0.31% Glutathione, 0.79% Tris HCl

## General Info

<b>Relevance</b>	PCQAP / MED15 is a member of the mediator complex that can act as both a positive and negative regulator of transcription. The mediator complex interacts with RNA pol II and its associated proteins, the general transcription factors. Mediator proteins serve as a bridge connecting the activator and basal transcription machinery.
------------------	---

<b>Cellular localization</b>	Cytoplasmic and Nuclear
------------------------------	-------------------------

## Images



ab162600 on a 12.5% SDS-PAGE stained with Coomassie Blue.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors