

Product datasheet

Recombinant Human PHIP protein ab196110

1 Image

Description

Product name	Recombinant Human PHIP protein
Purity	> 90 % SDS-PAGE.
Expression system	Escherichia coli
Accession	<u>Q8WWQ0</u>
Protein length	Protein fragment
Animal free	No
Nature	Recombinant
Species	Human
Sequence	MHHHHHHSLYKPLDGEWGTNPRDEECERIVAGINQLMTL DIASAFVAPV DLQAYPMYCTVVAYPTDLSTIKQRLNRFYRRVSSLMWEV RYIEHNTRTF NEPGSPIMKSAKFVTDLLLHFIKDQTCYNIPLYNMKKKKVL SDSEDEE
Predicted molecular weight	17 kDa including tags
Amino acids	1146 to 1287
Tags	His tag N-Terminus
Additional sequence information	Bromodomain 1 (NM_017934).

Specifications

Our **Abpromise guarantee** covers the use of **ab196110** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications SDS-PAGE

Form Liquid

Preparation and Storage

Stability and Storage Shipped on Dry Ice. Store at -80°C.
pH: 8.00


Constituents: 0.63% Tris HCl, 0.64% Sodium chloride, 0.02% Potassium chloride, 20% Glycerol (glycerin, glycerine)

General Info

Relevance

PHIP contains two bromo domains and eight WD repeats. It may be involved in insulin signaling. It is expressed in myeloma and epidermoid carcinoma cell lines.

Images



4-20% SDS-PAGE analysis of 2.3 µg ab196110 with Coomassie staining.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors