

Product datasheet

Recombinant Human RPL34 protein (denatured)  
ab171711

1 Image

Description	
Product name	Recombinant Human RPL34 protein (denatured)
Purity	> 85 % SDS-PAGE.
Expression system	Escherichia coli
Accession	<u>P49207</u>
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	MGSSHHHHHHSSGLVPRGSHMGSMVQRLTYRRRLSYNT ASNKTRLSRTPG NRIVLYTKKVGKAPKSACGVCPGRLRGVRAVRPKVLMR LSKTKKHVSRA YGGSMCAKCVRDRIKRAFLIEEQKIVVKVLKAQAQSQKAK
Predicted molecular weight	16 kDa including tags
Amino acids	1 to 117
Tags	His tag N-Terminus
Description	Recombinant Human RPL34 protein

Specifications	
Our <b>Abpromise guarantee</b> covers the use of <b>ab171711</b> in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
Applications	SDS-PAGE
Form	Liquid

Preparation and Storage	
Stability and Storage	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -

80°C. Avoid freeze / thaw cycle.

pH: 8.00

Constituents: 0.32% Tris HCl, 10% Glycerol (glycerin, glycerine), 2.34% Sodium chloride

## General Info

---


### Sequence similarities

Belongs to the ribosomal protein L34e family.

---

## Images

---



15% SDS-PAGE analysis of ab171711 (3µg).

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors