

Recombinant Human RPP30 protein ab102557

1 Image

Description	
Product name	Recombinant Human RPP30 protein
Purity	> 90 % SDS-PAGE. ab102557 was purified using conventional chromatography techniques.
Expression system	Escherichia coli
Accession	<b><u>P78346</u></b>
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	<b>MGSSHHHHHHSSGLVPRGSHMG</b> SMAVFADLDLRAGSD LKALRGLVETAAH LGYSVVAINHIVDFKEKKQEIEKPVAVSELFTTLPVQGKSR PIKILTRL TIIVSDPSHCNVL RATSSRARLYDVVAVFPKTEKLFHIACTH LDVDLVCI TVTEKLPHYFKRPPINVAIDRGLAFELVYSPAIKDSTMRRYTI SSALNLM QICKGKNVISSAAERPLEIRGPYDVANLGLLFGLSESDAKA AVSTNCRA ALLHGETRKTAFGIISTVKKPRPSEGDEDCLPASKKAKCE G
Predicted molecular weight	32 kDa including tags
Amino acids	1 to 268
Tags	His tag N-Terminus

Specifications	
Our <b><u>Abpromise guarantee</u></b> covers the use of <b>ab102557</b> in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
Applications	Mass Spectrometry SDS-PAGE
Mass spectrometry	MALDI-TOF

**Form** Liquid

## Preparation and Storage

### Stability and Storage


Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

pH: 8.00

Constituents: 0.077% DTT, 0.316% Tris HCl, 0.0292% EDTA, 20% Glycerol (glycerin, glycerine), 1.16% Sodium chloride

## General Info

## Images



15% SDS-PAGE analysis of ab102557 (3 µg).

SDS-PAGE - Recombinant Human RPP30 protein  
(ab102557)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors