

Recombinant Human RTR protein ab114813

1 Image

Description

Product name	Recombinant Human RTR protein		
Expression system	Wheat germ		
Accession	<u>Q15406</u>		
Protein length	Full length protein		
Animal free	No		
Nature	Recombinant		
Species	Human		
Sequence	MERDEPPPSGGGGGGGSAGFLEPPAALPPPPRNGFCQ DELAELDPGTNDR AEQRTCLICGDRATGLHYGIISCEGCKGFFKRSICNKRKYR CSRDKNCVM SRKQRNRCQYCRLLKCLQMGMNRKAIREDGMPGGRNKS IPVQISEEEIE RIMSGQEFEEEEANHWSNHGSDSHSSPGNRASESNQPSP GSTLSSRSVELN GFMAFREQYMGMSVPPHYQYIPHLFSYSGHSPLLPQAR SLDPQSYSLIH QLLSAEDLE PLGTPMLIED GYAVTQAE LFLCRLADEL LFRQIAWIKK LPFFCELSIK DYTCLLSSTW QELILLSSLT VYSKQIFGEL ADVTAKYSPS DEELHRFSDE GMEVIERLIY LYHKFHLKLV SNEEYACMKA INFLNQDIRG LTSASQLEQL NKRYWYICQD FTEYKYTHQP NRFPDLMMCL PEIRYAGKM VNVPLEQLPL LFKVVVLSCK TSVGKE		
Predicted molecular weight	81 kDa including tags		
Amino acids	1 to 475		

Specifications

Our **Abpromise guarantee** covers the use of **ab114813** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	ELISA
	SDS-PAGE

	Western blot
Form	Liquid


Preparation and Storage

Stability and Storage	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00 Constituents: 0.3% Glutathione, 0.79% Tris HCl
------------------------------	---

General Info

Function	Orphan nuclear receptor. Binds to a response element containing the sequence 5'-TCAAGGTCA-3'. May be involved in the regulation of gene expression in germ cell development during gametogenesis.
Tissue specificity	Shows highest expression in the germ cells of the adult testis.
Sequence similarities	Belongs to the nuclear hormone receptor family. NR6 subfamily. Contains 1 nuclear receptor DNA-binding domain.
Cellular localization	Nucleus.

Images



12.5% SDS-PAGE Stained with Coomassie Blue

SDS-PAGE - Recombinant Human RTR protein
(ab114813)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors