

Recombinant human TIMP1 protein ab82104

★★★★★ [2 Abreviews](#) [1 References](#) [2 Images](#)

Description

Product name	Recombinant human TIMP1 protein
Biological activity	Biological Activity : TIMP1 activity was measured by its ability to inhibit human MMP1 induced hydrolysis of a chromogenic peptide substrate at room temperature. Half maximal inhibition was obtained at a TIMP1 concentration of approximately 0.5 µg/ml, when using an MMP1 concentration of 1.6 µg/ml.
Purity	> 95 % SDS-PAGE. Purity is greater than 95% by SDS-PAGE gel and HPLC analyses. Endotoxin level is less than 0.1 ng per µg (1EU/µg).
Expression system	Escherichia coli
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	CTCVPPHPQT AFCNSDLVIR AKFVGTPPEVN QTTLYQRYEI KMTKMYKGFQ ALGDAADIRF VYTPAMESVC GYFHRSHNRS EEFLIAGKLQ DGLLHITTCF FVAPWNSLSL AQRGFTKTY TVGCEECTVF PCLSIPCKLQ SGTHCLWTDQ LLQGSEKGFQ SRHLACLPRE PGLCTWQSLR SQIA

Specifications

Our **Abpromise guarantee** covers the use of **ab82104** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	Functional Studies SDS-PAGE Western blot
Form	Lyophilized
Additional notes	Centrifuge the vial prior to opening. Reconstitute in water to a concentration of 0.1-1.0 mg/ml. Do not vortex (can be stored at 4C for up to 1 week). For extended storage, it is recommended to further dilute in a buffer containing a carrier protein

(example 0.1% BSA) and store in working aliquots at -20°C to -80°C.

Preparation and Storage

Stability and Storage

Shipped at 4°C. The lyophilized protein is stable for a few weeks at room temperature. Store at -20°C or -80°C. Avoid freeze / thaw cycle. Please see notes section.

This product is an active protein and may elicit a biological response in vivo, handle with caution.

Reconstitution

Centrifuge the vial prior to opening. Reconstitute in water to a concentration of 0.1-1.0 mg/ml. Do not vortex (can be stored at 4°C for up to 1 week). For extended storage, it is recommended to further dilute in a buffer containing a carrier protein (example 0.1% BSA) and store in working aliquots at -20°C to -80°C.

General Info

Function

Complexes with metalloproteinases (such as collagenases) and irreversibly inactivates them by binding to their catalytic zinc cofactor. Also mediates erythropoiesis in vitro; but, unlike IL-3, it is species-specific, stimulating the growth and differentiation of only human and murine erythroid progenitors. Known to act on MMP-1, MMP-2, MMP-3, MMP-7, MMP-8, MMP-9, MMP-10, MMP-11, MMP-12, MMP-13 and MMP-16. Does not act on MMP-14.

Sequence similarities

Belongs to the protease inhibitor I35 (TIMP) family.

Contains 1 NTR domain.


Post-translational modifications

The activity of TIMP1 is dependent on the presence of disulfide bonds.

Cellular localization

Secreted.

Images



Sandwich ELISA - Recombinant human TIMP1 protein (ab82104)



Western blot - Recombinant human TIMP1 protein
(ab82104)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors