

Recombinant Human Transmembrane protein 30A
ab163206

1 Image

Description

Product name	Recombinant Human Transmembrane protein 30A		
Expression system	Wheat germ		
Accession	<u>Q9NV96</u>		
Protein length	Full length protein		
Animal free	No		
Nature	Recombinant		
Species	Human		
Sequence	MAMNYNAKDEV DGGPPCAPGGTAKTRPDNTAFKQQRL PAWQPILTAGTV LPIFFIIGLIFIPIGIGIFVTSNNIREIEIDYTGTEPSSPCNKCLSP DVT PCFCTINF TLEKSFEGNVFMYG LSNFYQNHRRYVKS RDD SQLNGDSSAL LNPSKECEPYRRNEDKPIAPCGA IANSMFNDTLELFLIGND SYIPIALK KKGIAWWTDKNVKFRNPPGGDNLEERFKGTTKPVNWLK PVYMLDSDPDNN GFINEDFIVWMRTAALPTFRKLYRLIERKSDLHPTLPAGRYS LNVTYNYP VHYFDGRKRMILSTISWMGGKNPFLGIAYI AVGSISFLLGVV LLVINHKYRNSSNTADITI		
Amino acids	1 to 361		
Tags	GST tag N-Terminus		

Specifications

Our Abpromise guarantee covers the use of ab163206 in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
Applications	ELISA Western blot
Form	Liquid

Additional notes

Preparation and Storage

Stability and Storage

Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.

pH: 8.00

Constituents: 0.31% Glutathione, 0.79% Tris HCl

General Info


Relevance

Belongs to the CDC50/LEM3 family.

Cellular localization

Membrane; Multi-pass membrane protein

Images



ab163206 on a 12.5% SDS-PAGE stained with Coomassie Blue.

SDS-PAGE - Recombinant Human Transmembrane protein 30A (ab163206)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors