

Recombinant Human TSC-1 protein ab54406

Description

Product name	Recombinant Human TSC-1 protein
Purity	> 95 % SDS-PAGE. > 98% by SDS-PAGE and HPLC analyses Endotoxin level: < 0.1 ng per ug of TSC-1
Expression system	Escherichia coli
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human

Specifications

Our **Abpromise guarantee** covers the use of **ab54406** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	SDS-PAGE Functional Studies
Form	Lyophilized
Additional notes	Previously labelled as Eotaxin 3.

Preparation and Storage

Stability and Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. Constituent: 0.121% Tris
Reconstitution	This protein is soluble in water and most aqueous buffers. The lyophilized protein can be reconstituted in water to a concentration of 0.1-1.0 mg/ml. It is recommended that further dilutions be made into buffer containing carrier protein or medium containing serum.

General Info

Function	Chemotactic for eosinophils and basophils. Binds to CCR3.
Tissue specificity	Ubiquitously expressed at low levels in various tissues including heart and ovary.

Sequence similarities	Belongs to the intercrine beta (chemokine CC) family.
Cellular localization	Secreted.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors