

Product datasheet

Recombinant Human TSR2 protein ab183259

1 Image

Description

Product name	Recombinant Human TSR2 protein
Purity	> 90 % SDS-PAGE. ab183259 was purified by using anion-exchange chromatography (DEAE sepharose resin) and gel-filtration chromatography (Sephacryl S-200) with 20mM Tris pH 7.5, 2mM EDTA.
Expression system	Escherichia coli
Accession	<u>Q969E8</u>
Protein length	Protein fragment
Animal free	No
Nature	Recombinant
Species	Human
Sequence	MGSSHHHHHSSGLVPRGSHMGSMAGAAEDARALFRA GVCAALEAWPALQ IAVENFGGVHSQEKAKWLGGAVEDYFMRNADLELDEV EDFLGELLTNEF DTVVEDGSLPQVSQQLQTMFHHFQRGDGAALREMASCIT QRKCKVTATAL KTARETDEDDVDSVEEMEVTATNDGAATDGVCPQPEP SDPDAQTIKEED IVEDGWTIVRRK
Predicted molecular weight	23 kDa including tags
Amino acids	1 to 190
Tags	His tag N-Terminus
Additional sequence information	NP_477511

Specifications

Our **Abpromise guarantee** covers the use of **ab183259** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	Mass Spectrometry SDS-PAGE
Mass spectrometry	MALDI-TOF
Form	Liquid

Preparation and Storage

Stability and Storage


Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 7.40

Constituents: 90% PBS, 10% Glycerol (glycerin, glycerine)

General Info

Images



15% SDS-PAGE analysis of ab183259, 3µg loaded. N.B. 23.2kDa (confirmed by MALDI-TOF), but molecular size on SDS-PAGE appears larger.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors