

Product datasheet

Recombinant Human TSTD3 protein ab167916

1 Image

Description

| | |
|-----------------------------------|--|
| Product name | Recombinant Human TSTD3 protein |
| Purity | > 90 % SDS-PAGE. purified by using conventional chromatography techniques. |
| Expression system | Escherichia coli |
| Accession | <u>H0UI37</u> |
| Protein length | Full length protein |
| Animal free | No |
| Nature | Recombinant |
| Species | Human |
| Sequence | MGSSHHHHHH SSGLVPRGSH MGSMKIEKCG WSEGLTSIKG NCHNFYTAIS KDVTYKELKN LLNSKNIMLI DVREIWEILE YQKIPESINV PLDEVGEALQ MNPRDFKEKY NEVKPSKSDS |
| Predicted molecular weight | 14 kDa including tags |
| Amino acids | 1 to 97 |
| Tags | His tag N-Terminus |

Specifications

Our **Abpromise guarantee** covers the use of **ab167916** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| | |
|--------------------------|-------------------------------|
| Applications | Mass Spectrometry SDS-PAGE |
| Mass spectrometry | MALDI-TOF |
| Form | Liquid |

Preparation and Storage

| | |
|------------------------------|--|
| Stability and Storage | Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle. |
|------------------------------|--|

pH: 8.00


Constituents: 0.02% DTT, 0.32% Tris HCl, 10% Glycerol (glycerin, glycerine), 0.88% Sodium chloride

General Info

Sequence similarities

Contains 1 rhodanese domain.

Images



15% SDS-PAGE analysis of 3 µg ab167916

SDS-PAGE - Recombinant Human TSTD3 protein
(ab167916)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors