

Product datasheet

Recombinant Human Utrophin protein ab159789

1 Image

Description

Product name	Recombinant Human Utrophin protein
Expression system	Wheat germ
Protein length	Protein fragment
Animal free	No
Nature	Recombinant
Species	Human
Sequence	LEARMQILEDHKNQLESQQLHRLRQLLEQPESDSRINGVSP WASPQHSALS YSLDPDASGPQFHQAAGEDLLAPPHDTSTDLTEVMEQIH STFPSCCPNVP SRPQAM
Amino acids	3328 to 3433
Tags	GST tag N-Terminus

Specifications

Our **Abpromise guarantee** covers the use of **ab159789** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	Western blot ELISA
Form	Liquid
Additional notes	

Preparation and Storage


Stability and Storage	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00 Constituents: 0.31% Glutathione, 0.79% Tris HCl
------------------------------	--

General Info

General info

Function	May play a role in anchoring the cytoskeleton to the plasma membrane.
Tissue specificity	Muscle.
Sequence similarities	Contains 2 CH (calponin-homology) domains. Contains 20 spectrin repeats. Contains 1 WW domain. Contains 1 ZZ-type zinc finger.
Cellular localization	Cell junction > synapse > postsynaptic cell membrane. Cytoplasm > cytoskeleton. Neuromuscular junction.

Images



ab159789 on a 12.5% SDS-PAGE stained with Coomassie Blue.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors