

# Recombinant Rubella Virus gpE1 protein ab43032

### Description

<b>Product name</b>	Recombinant Rubella Virus gpE1 protein
<b>Purity</b>	> 95 % SDS-PAGE. Purified by ion exchange chromatography.
<b>Expression system</b>	Escherichia coli
<b>Protein length</b>	Protein fragment
<b>Animal free</b>	No
<b>Nature</b>	Recombinant

### Specifications

Our **Abpromise guarantee** covers the use of **ab43032** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

<b>Applications</b>	ELISA Western blot
<b>Form</b>	Liquid

### Preparation and Storage

<b>Stability and Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. Preservative: 0.14% Imidazole Constituents: 48% Urea, 1.74% Sodium chloride
------------------------------	---

### General Info

<b>Relevance</b>	Rubella virus is the only member of the Rubrivirus genus of the Togavirus family. Unlike most Togaviruses it is NOT arthropod borne, but is acquired via the respiratory route. It is an enveloped (toga=cloak), non-segmented, positive sense, RNA virus and replicates in the cytoplasm. It consists of 3 structural proteins; E1,E2 membrane bound glycoproteins, and C capsid protein. Rubella or German measles is a mild disease in children and adults, but can cause devastating problems if it infects the fetus, especially if infection is in the first few weeks of pregnancy.
------------------	--

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors