

Recombinant VP16 protein ab81878

Description

Product name	Recombinant VP16 protein
Purity	> 95 % SDS-PAGE. ab81878 is isolated from an E. coli strain that carries the coding sequence of the fused protein under the control of a T7 promoter.
Expression system	Escherichia coli
Protein length	Protein fragment
Animal free	No
Nature	Recombinant
Amino acids	1 to 147

Specifications

Our **Abpromise guarantee** covers the use of **ab81878** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	SDS-PAGE
Form	Liquid

Preparation and Storage

Stability and Storage	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 7.9 Constituents: 0.75% Potassium chloride, 0.0154% DTT, 0.316% Tris HCl, 0.00584% EDTA, 20% Glycerol (glycerin, glycerine)
------------------------------	---

General Info

Relevance	VP16 belongs to the herpesviridae alpha trans-inducing protein family and is responsible for transcriptional activation of immediate early promoters (alpha genes). VP16 has been shown to bind to TBP, TFIIIB, and replication factor A.
------------------	---

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors