



Product datasheet

25x PBS Buffer pH 7.6 ab64026

1 References

Overview

Product name	25x PBS Buffer pH 7.6
Tested applications	Suitable for: IHC-P
General notes	In order to provide this product at 25X, the salts in solution are at a very high concentration. It has been found that at temperatures below 44 degrees F (6 degrees C), the salts fall out of solution. If this occurs during your shipment conditions, please stir the product to get the salts back into solution:

Properties

Form	Liquid
Storage instructions	Store at room temperature.
Storage buffer	pH: 7.60 Constituent: PBS Buffer 25x concentrated.
Relevance	PBS (Phosphate Buffered Saline) is used with immunolabeling techniques as a wash buffer.

Applications

The Abpromise guarantee Our **Abpromise guarantee** covers the use of ab64026 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
IHC-P		Use at an assay dependent dilution.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors