

## Product datasheet

# Donkey Anti-Goat IgG H&L (HRP) preadsorbed ab7125

[14 References](#) [1 Image](#)

### Overview

---

<b>Product name</b>	Donkey Anti-Goat IgG H&L (HRP) preadsorbed	
<b>Host species</b>	Donkey	
<b>Target species</b>	Goat	
<b>Tested applications</b>	<b>Suitable for:</b> Dot blot, ELISA, IHC-P, IHC-Fr, Immunomicroscopy, WB, ICC/IF	
<b>Minimal cross-reactivity</b>	Chicken, Guinea pig, Hamster, Horse, Mouse, Rabbit, Rat	<a href="#">more details</a>
<b>Immunogen</b>	Full length native Goat IgG (purified).	
<b>Conjugation</b>	HRP	

### Properties

---

<b>Form</b>	Liquid
<b>Storage instructions</b>	Shipped at 4°C. Store at +4°C.
<b>Storage buffer</b>	pH: 6.50 Preservative: 0.01% Gentamicin sulphate Constituents: 0.42% Potassium phosphate, 0.87% Sodium chloride, 1% BSA
<b>Purity</b>	Affinity purified
<b>Purification notes</b>	This product was prepared from monospecific antiserum by immunoaffinity chromatography using Goat IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities.
<b>Conjugation notes</b>	Horseradish Peroxidase (HRP)
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG


### Applications

---

**The Abpromise guarantee** Our [Abpromise guarantee](#) covers the use of ab7125 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
Dot blot		Use at an assay dependent dilution.
ELISA		1/170000.
IHC-P		Use at an assay dependent dilution.
IHC-Fr		Use at an assay dependent dilution.
Immunomicroscopy		Use at an assay dependent dilution.
WB		1/2000 - 1/10000.
ICC/IF		1/500 - 1/2500.

## Images



ab7125 was used at dilution 1/10000 with the primary antibody [ab37917](#) in WB. See the review on [ab37917](#).

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise,

please visit <https://www.abcam.com/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors