

Product datasheet

Donkey Anti-Rabbit IgG H&L (12nm Gold) preadsorbed ab105295

1 References

Overview

Product name	Donkey Anti-Rabbit IgG H&L (12nm Gold) preadsorbed
Host species	Donkey
Target species	Rabbit
Specificity	Recognizes rabbit IgG, H&L. X-Adsorbed and exhibits minimal crossreactivity to bovine, chicken, goat, guinea pig, Syrian hamster, equine, human, mouse, rat and sheep serum proteins.
Tested applications	Suitable for: Electron Microscopy, ELISA, WB
Minimal cross-reactivity	Chicken, Cow, Goat, Guinea pig, Horse, Human, Mouse, Rat, Sheep, Syrian hamster
Conjugation	Gold 12nm

[more details](#)

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C.
Storage buffer	pH: 8.50 Preservative: 0.05% Sodium azide Constituents: 0.164% Sodium phosphate, 0.381% Sodium borate, 1.5% BSA, 0.87% Sodium chloride
Purity	Immunogen affinity purified
Clonality	Polyclonal
Isotype	IgG

Applications

The Abpromise guarantee Our [Abpromise guarantee](#) covers the use of ab105295 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
Electron Microscopy		1/20 - 1/40.
ELISA		Use at an assay dependent dilution.
WB		Use at an assay dependent dilution.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors