

Product datasheet

Donkey Anti-Rabbit IgG H&L (6nm Gold) preadsorbed ab105294

1 References

Overview

Product name	Donkey Anti-Rabbit IgG H&L (6nm Gold) preadsorbed
Host species	Donkey
Target species	Rabbit
Specificity	Recognizes the heavy chains on rabbit IgG and the light chains common to most rabbit immunoglobulins. Does not crossreact with non-immunoglobulin serum proteins. Exhibits minimal crossreactivity to bovine, chicken, goat, guinea pig, syrian hamster, horse, human, mouse, rat and sheep serum proteins.
Tested applications	Suitable for: Electron Microscopy, WB, ELISA
Minimal cross-reactivity	Chicken, Cow, Goat, Guinea pig, Horse, Human, Mouse, Rat, Sheep, Syrian hamster
Conjugation	Gold 6nm

[more details](#)

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C.
Storage buffer	pH: 8.50 Preservative: 0.05% Sodium azide Constituents: 0.164% Sodium phosphate, 0.381% Sodium borate, 1.5% BSA, 0.87% Sodium chloride
Purity	Immunogen affinity purified
Clonality	Polyclonal
Isotype	IgG

Applications

The Abpromise guarantee Our **Abpromise guarantee** covers the use of ab105294 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
Electron Microscopy		1/20 - 1/40.
WB		Use at an assay dependent concentration.
ELISA		Use at an assay dependent concentration.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors