

Product datasheet

Donkey Anti-Rabbit IgG H&L (Texas Red ®) preadsorbed ab7081

5 References

Overview

Product name	Donkey Anti-Rabbit IgG H&L (Texas Red ®) preadsorbed
Host species	Donkey
Target species	Rabbit
Tested applications	Suitable for: ICC/IF, Immunomicroscopy, Flow Cyt, IHC-P, IHC-Fr, ELISA
Minimal cross-reactivity	Chicken, Cow, Goat, Guinea pig, Hamster, Horse, Human, Mouse, Rat, Sheep more details
Immunogen	Full length native Rabbit IgG (purified).
Conjugation	Texas Red ®. Ex: 596nm, Em: 620nm

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
Storage buffer	pH: 6.50 Preservative: 0.01% Sodium azide Constituents: 0.42% Tripotassium orthophosphate, 0.87% Sodium chloride, 1% BSA
Purity	Affinity purified
Purification notes	Texas Red Sulfonyl Chloride (Molecular Weight 625 daltons) Absorption Wavelength: 596 nm Emission Wavelength: 620 nm Fluorochrome/Protein Ratio: 1.7 moles Texas Red per mole of Goat IgG
Conjugation notes	This product was prepared from monospecific antiserum by immunoaffinity chromatography using Rabbit IgG coupled to agarose
Clonality	Polyclonal

Applications

The **Abpromise guarantee** Our [Abpromise guarantee](#) covers the use of ab7081 in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
ICC/IF		
Immunomicroscopy		
Flow Cyt		
IHC-P		
IHC-Fr		
ELISA		

Application notes

Flow Cyt: Use at an assay dependant concentration.
IHC-P: Use at an assay dependant concentration.
IHC-Fr: Use at an assay dependant concentration.
IM: Use at an assay dependant concentration.

Not tested in other applications.
Optimal dilutions/concentrations should be determined by the end user.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors