

## Product datasheet

### Goat Anti-Chicken IgY H&L (Alkaline Phosphatase) preadsorbed ab7119

#### Overview

<b>Product name</b>	Goat Anti-Chicken IgY H&L (Alkaline Phosphatase) preadsorbed	
<b>Host species</b>	Goat	
<b>Target species</b>	Chicken	
<b>Tested applications</b>	<b>Suitable for:</b> ICC/IF, Dot blot, ELISA, IHC-P, IHC-Fr, Immunomicroscopy, WB	
<b>Minimal cross-reactivity</b>	Cow, Goat, Guinea pig, Hamster, Horse, Human, Rabbit, Rat, Sheep	<a href="#">more details</a>
<b>Immunogen</b>	Full length native chicken IgY (purified).	
<b>Conjugation</b>	Alkaline Phosphatase	

#### Properties

<b>Form</b>	Liquid
<b>Storage instructions</b>	Shipped at 4°C. Store at +4°C.
<b>Storage buffer</b>	Constituents: 0.00136% Zinc chloride, 0.0095% Magnesium chloride, 0.79% Tris HCl, 50% Glycerol (glycerin, glycerine), 0.87% Sodium chloride
<b>Purity</b>	Affinity purified
<b>Purification notes</b>	This product was prepared from monospecific antiserum by immunoaffinity chromatography using chicken IgG coupled to agarose beads.
<b>Conjugation notes</b>	Alkaline Phosphatase (Calf Intestine) (Molecular Weight 140,000 daltons)
<b>Clonality</b>	Polyclonal

#### Applications

**The Abpromise guarantee** Our [Abpromise guarantee](#) covers the use of ab7119 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
ICC/IF		

Application	Abreviews	Notes
Dot blot		
ELISA		
IHC-P		
IHC-Fr		
Immunomicroscopy		
WB		

#### Application notes

Dot: Use at an assay dependant dilution.

ELISA: 1/1000 - 1/4000.

IHC-P: Use at an assay dependant dilution.

IHC-Fr: Use at an assay dependant dilution.

IM: Use at an assay dependant dilution.

WB: Use at an assay dependant dilution.

Not tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

#### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

#### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors