

Product datasheet

Goat F(ab')₂ Anti-Mouse IgG - H&L (DyLight® 488), pre-adsorbed ab98726

★★★★★ [1 Abreviews](#) [1 References](#) [1 Image](#)

Overview

Product name	Goat F(ab') ₂ Anti-Mouse IgG - H&L (DyLight® 488), pre-adsorbed
Host species	Goat
Target species	Mouse
Tested applications	Suitable for: IHC-P, ICC/IF, Flow Cyt
Minimal cross-reactivity	Human, Rat more details
Immunogen	Other Immunogen Type corresponding to IgG.
Conjugation	DyLight® 488. Ex: 493nm, Em: 518nm

Properties


Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C.
Storage buffer	pH: 6.8 Preservative: 0.09% Sodium azide Constituents: PBS, 0.2% BSA
Purity	Immunogen affinity purified
Purification notes	This antibody was isolated by affinity chromatography using antigen coupled to agarose beads. F(ab') ₂ fragment were generated using a pepsin digestion. Fc fragments and whole IgG molecules have been removed. Fragments were conjugated to DyLight® 488.
Clonality	Polyclonal
Isotype	IgG

Applications

The Abpromise guarantee Our **Abpromise guarantee** covers the use of ab98726 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
IHC-P		1/50 - 1/500.
ICC/IF		1/50 - 1/500.
Flow Cyt	★★★★★ (1)	1/50 - 1/200.

Images



Emission spectra of DyLight® fluorochromes available in our catalog.

Line colors represent the approximate visible colors of the wavelength maxima.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors