

Product datasheet

Goat Anti-Human IgG+IgM+IgA H&L (Biotin) ab102418

1 References

Overview

| | |
|----------------------------|---|
| Product name | Goat Anti-Human IgG+IgM+IgA H&L (Biotin) |
| Host species | Goat |
| Target species | Human |
| Specificity | By immunoelectrophoresis and ELISA this antibody reacts with human IgG, IgA and IgM. No antibody was detected against non-immunoglobulin serum proteins. This antibody may cross react with immunoglobulins from other species. |
| Tested applications | Suitable for: IHC-P, ICC/IF, ELISA, WB |
| Immunogen | Other Immunogen Type corresponding to IgG + IgM + IgA. |
| Conjugation | Biotin |

Properties

| | |
|-----------------------------|--|
| Form | Liquid |
| Storage instructions | Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle. |
| Storage buffer | pH: 6.8 Preservative: 0.09% Sodium azide Constituents: PBS, 0.2% BSA |
| Purity | Immunogen affinity purified |
| Purification notes | This antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to Biotin. |
| Clonality | Polyclonal |
| Isotype | IgG |

Applications

The Abpromise guarantee Our **Abpromise guarantee** covers the use of ab102418 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| Application | Abreviews | Notes |
|-------------|-----------|---|
| IHC-P | | 1/500 - 1/5000. |
| ICC/IF | | 1/200 - 1/1000. |
| ELISA | | 1/10000 - 1/200000. ELISA - Primary |
| WB | | 1/10000 - 1/200000. WB - Chemiluminescent |

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors