

Product datasheet

Goat Anti-Human IgM mu chain (PE) preadsorbed
ab178669

1 Image

Overview

Product name	Goat Anti-Human IgM mu chain (PE) preadsorbed	
Host species	Goat	
Target species	Human	
Tested applications	Suitable for: Flow Cyt	
Minimal cross-reactivity	Cow	more details
Conjugation	PE. Ex: 488nm, Em: 575nm	

Properties


Form	Lyophilized:Reconstitute with 1 ml of distilled water
Storage instructions	Shipped at 4°C. Store at +4°C. Do Not Freeze. Store In the Dark.
Storage buffer	pH: 7.40 Preservative: 0.05% Sodium azide Constituents: 1.4% Sodium phosphate, 0.58% Sodium chloride, 0.1% BSA, 1.71% Sucrose
Clonality	Polyclonal
Isotype	IgG
General notes	Product should be stored at 2-8 degC in the dark and be used within 6 months. If further dilution of the conjugate is required, use diluted material within one week.

Applications

The Abpromise guarantee Our **Abpromise guarantee** covers the use of ab178669 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
Flow Cyt		Use at an assay dependent concentration.

Images



Fluorescence excitation and emission spectra of R-phycoerythrin in 100 mM sodium phosphate (pH 7.2) + 1 mM EDTA and 1 mM sodium azide. Emission scan was taken with excitation at 498 nm. Excitation scan was taken with emission at 575 nm. Scans were normalized to equalize peak heights.

Flow Cytometry - Goat Anti-Human IgM mu chain
(Phycoerythrin) preadsorbed (ab178669)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors