



Product datasheet

Goat Anti-Mouse IgG1 heavy chain (FITC) preadsorbed ab98692

3 References

Overview

Product name	Goat Anti-Mouse IgG1 heavy chain (FITC) preadsorbed		
Host species	Goat		
Target species	Mouse		
Tested applications	Suitable for: IHC-P, ICC/IF, Flow Cyt		
Minimal cross-reactivity	Human, Rat		more details
Conjugation	FITC. Ex: 493nm, Em: 528nm		

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C.
Storage buffer	pH: 6.8 Preservative: 0.09% Sodium azide Constituents: PBS, 0.2% BSA
Purity	Immunogen affinity purified
Purification notes	ab98692 was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to FITC.
Clonality	Polyclonal
Isotype	IgG

Applications

The Abpromise guarantee Our **Abpromise guarantee** covers the use of ab98692 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
IHC-P		1/50 - 1/500.

Application	Abreviews	Notes
ICC/IF		1/50 - 1/500.
Flow Cyt		1/50 - 1/200.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors