

Goat Anti-Mouse IgG1 H&L (PE) ab99912

Overview

Product name	Goat Anti-Mouse IgG1 H&L (PE)
Host species	Goat
Target species	Mouse
Tested applications	Suitable for: Flow Cyt
Conjugation	PE. Ex: 488nm, Em: 575nm

Properties

Form	Lyophilized:rehydrate with 1.0 ml distilled deionized water
Storage instructions	Shipped at 4°C. Store at +4°C. Store In the Dark.
Storage buffer	pH: 7.40 Preservative: 0.05% Sodium azide Constituents: 6.84% Sucrose, 0.1% BSA, 0.87% Sodium chloride, 1.64% Sodium phosphate
Purity	Protein G purified
Conjugation notes	R-phycoerythrin Excitation max.: 565>498 nm Emission max.: 578 nm
Clonality	Polyclonal
Isotype	IgG

Applications

The Abpromise guarantee Our **Abpromise guarantee** covers the use of ab99912 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
Flow Cyt		Use at an assay dependent concentration.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors