

Product datasheet

Goat Anti-Mouse IgG+IgM H&L (HRP) preadsorbed ab47827

[26 References](#) [1 Image](#)

Overview

| | |
|---------------------------------|---|
| Product name | Goat Anti-Mouse IgG+IgM H&L (HRP) preadsorbed |
| Host species | Goat |
| Target species | Mouse |
| Specificity | Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Goat Serum, anti-Peroxidase, Mouse IgG and Mouse IgM. |
| Tested applications | Suitable for: ICC/IF, ELISA, WB, IHC-P, IHC-Fr, Electron Microscopy |
| Minimal cross-reactivity | Cow, Horse, Human more details |
| Immunogen | Full length native proteins (purified) (Mouse) |
| Conjugation | HRP |

Properties

| | |
|-----------------------------|---|
| Form | Liquid |
| Storage instructions | Shipped at 4°C. Store at +4°C. |
| Storage buffer | pH: 6.50 Preservative: 0.01% Gentamicin sulphate Constituents: 0.42% Tripotassium orthophosphate, 0.87% Sodium chloride, 1% BSA |
| Purity | IgG and protease free Immunogen affinity purified |
| Purification notes | Prepared from polyspecific antiserum by immunoaffinity chromatography using antigens coupled to agarose beads followed by solid phase adsorption (s) to remove any unwanted reactivities. |
| Clonality | Polyclonal |
| Isotype | IgG |

Applications


The Abpromise guarantee

Our **Abpromise guarantee** covers the use of ab47827 in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| Application | Abreviews | Notes |
|---------------------|-----------|-------------------------------------|
| ICC/IF | | 1/500 - 1/2500. |
| ELISA | | 1/5000 - 1/40000. |
| WB | | 1/2000 - 1/10000. |
| IHC-P | | 1/500 - 1/2500. |
| IHC-Fr | | 1/500 - 1/2500. |
| Electron Microscopy | | Use at an assay dependent dilution. |

Images



Ab47827 was used at dilution 1/500 with the primary antibody **ab2869** in IHC-P. See the review on **ab2869**.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Goat Anti-Mouse IgG+IgM H&L (HRP) preadsorbed (ab47827)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors