

Product datasheet

Goat Anti-Mouse IgM (SureLight® APC) ab99905

Overview

Product name	Goat Anti-Mouse IgM (SureLight® APC)
Host species	Goat
Target species	Mouse
Tested applications	Suitable for: Flow Cyt
Conjugation	SureLight® APC. Ex: 652nm, Em: 657nm

Properties

Form	Lyophilized: Lyophilized material is stable for up to one year. After product has been reconstituted, product should be stored at 2-8°C in the dark and be used within 3 months. If further dilution of the conjugate is required, use diluted material within one week.
Storage instructions	Shipped at 4°C. Store at +4°C.
Storage buffer	pH: 8.20 Constituents: 0.1% BSA, 0.87% Sodium chloride, 1.71% Sucrose
Purity	Protein G purified
Clonality	Polyclonal
Isotype	IgG

Applications

The Abpromise guarantee Our **Abpromise guarantee** covers the use of ab99905 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
Flow Cyt		Use a concentration of 3 µg/ml.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors