

Product datasheet

Goat anti-Rabbit IgG (SureLight® APC) ab72567

2 References

Overview

Product name	Goat anti-Rabbit IgG (SureLight® APC)
Host species	Goat
Target species	Rabbit
Tested applications	Suitable for: Flow Cyt, FRET
Minimal cross-reactivity	Cow, Mouse more details
Immunogen	Rabbit IgG
Conjugation	SureLight® APC. Ex: 652nm, Em: 657nm

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C.
Storage buffer	pH: 7.40 Preservative: 0.013% Sodium azide Constituents: 0.1% BSA, 0.87% Sodium chloride, 1.64% Sodium phosphate, 1.71% Sucrose
Purity	Size exclusion
Purification notes	Size exclusion chromatography as well as some proprietary methods for purifying the final materials were used.
Conjugation notes	Excitation max. 652 nm. Emission max. 657.5 nm. F/P 1:1.
Clonality	Polyclonal
Isotype	IgG

Applications

The Abpromise guarantee Our [Abpromise guarantee](#) covers the use of ab72567 in the following tested applications.
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
Flow Cyt		Use at an assay dependent concentration.
FRET		Use at an assay dependent concentration.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors