

Product datasheet

Goat Anti-Rat IgG H&L (Cy3 ®) preadsorbed ab98416

★★★★★ [1 Abreviews](#) [8 References](#)

Overview

Product name	Goat Anti-Rat IgG H&L (Cy3 ®) preadsorbed	
Host species	Goat	
Target species	Rat	
Tested applications	Suitable for: IHC-P, ICC/IF, Flow Cyt	
Minimal cross-reactivity	Chicken, Cow, Goat, Horse, Human, Mouse, Rabbit, Sheep	more details
Conjugation	Cy3 ®. Ex: 552nm, Em: 565nm	

Properties

Form	Liquid	
Storage instructions	Shipped at 4°C. Store at +4°C.	
Storage buffer	pH: 6.8	
	Preservative: 0.09% Sodium azide	
	Constituents: PBS, 0.2% BSA	
Purity	Immunogen affinity purified	
Purification notes	This antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to Cy3.	
Clonality	Polyclonal	
Isotype	IgG	

Applications

The Abpromise guarantee Our **Abpromise guarantee** covers the use of ab98416 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
IHC-P		1/50 - 1/500.

Application	Abreviews	Notes
ICC/IF		1/50 - 1/500.
Flow Cyt		1/100 - 1/1000.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors