

Anti-Human IgG4 pFc' [HP6023] ab99814

Overview

Product name	Anti-Human IgG4 pFc' [HP6023]
Host species	Mouse
Target species	Human
Specificity	ab99814 reacts with the pFc' portion of the heavy chain of Human IgG4.
Tested applications	Suitable for: ELISA, WB, Dot blot, IHC-Fr, ICC/IF
Immunogen	Full length protein corresponding to IgG4.

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
Storage buffer	pH: 8.20 Constituent: 100% Borate buffered saline No preservatives or amine-containing buffer salts added.
Purity	Protein A purified
Clonality	Monoclonal
Clone number	HP6023
Isotype	IgG3
Light chain type	kappa

Applications

The Abpromise guarantee Our **Abpromise guarantee** covers the use of ab99814 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
ELISA		Use at an assay dependent dilution.

Application	Abreviews	Notes
WB		Use at an assay dependent dilution.
Dot blot		Use at an assay dependent dilution. (for Dot and slot immunoblotting)
IHC-Fr		Use at an assay dependent dilution.
ICC/IF		Use at an assay dependent dilution.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors