

Product datasheet

Mouse monoclonal [JDC-10] Anti-Human IgG Fc ab99757

7 References

Overview

Product name	Mouse monoclonal [JDC-10] Anti-Human IgG Fc
Host species	Mouse
Target species	Human
Specificity	ab99757 reacts with the Fc portion of the heavy chain of all subclasses of Human IgG as demonstrated by ELISA; may also react with IgG from other species.
Tested applications	Suitable for: ICC/IF, ELISA, WB, Dot blot, ICC, IHC-Fr

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
Storage buffer	pH: 8.20 Constituent: 100% Borate buffered saline
	No preservatives or amine-containing buffer salts added.
Purity	Protein G purified
Clonality	Monoclonal
Clone number	JDC-10
Isotype	IgG1
Light chain type	kappa

Applications

The Abpromise guarantee Our **Abpromise guarantee** covers the use of ab99757 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
ICC/IF		Use 1µg for 10 ⁶ cells. for indirect immunofluorescent staining of IgG+ Human B lymphocytes.
ELISA		Use at an assay dependent dilution.
WB		Use at an assay dependent dilution.
Dot blot		Use at an assay dependent dilution.
ICC		Use at an assay dependent dilution.
IHC-Fr		Use at an assay dependent dilution.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors