

## Product datasheet

# Mouse monoclonal [MARM-4] Anti-Rat IgM H&L (HRP) ab11678

### Overview

Product name	Mouse monoclonal [MARM-4] Anti-Rat IgM H&L (HRP)
Host species	Mouse
Target species	Rat
Specificity	Ab11678 recognises rat IgM and does not react with other rat immunoglobulin isotypes.
Tested applications	<b>Suitable for:</b> ELISA
Immunogen	Full length protein corresponding to Rat IgM.
Conjugation	HRP

### Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C.
Storage buffer	Preservative: 0.09% Sodium azide Constituents: PBS, 50% Glycerol
Purification notes	Purified IgG prepared by affinity chromatography.
Clonality	Monoclonal
Clone number	MARM-4
Myeloma	Sp2/0
Isotype	IgG1

### Applications

**The Abpromise guarantee** Our **Abpromise guarantee** covers the use of ab11678 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
ELISA		

## Application notes

ELISA: Use at a concentration of 0.25 µg/ml as detecting antibody.

Not tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors