

Rabbit Anti-Chicken IgY H&L (Texas Red ®) ab6751

11 References

Overview

Product name	Rabbit Anti-Chicken IgY H&L (Texas Red ®)
Host species	Rabbit
Target species	Chicken
Tested applications	Suitable for: Immunomicroscopy, Flow Cyt, IHC-P, IHC-Fr, ICC/IF, ELISA
Immunogen	chicken IgY whole molecule
Conjugation	Texas Red ®. Ex: 596nm, Em: 620nm

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
Storage buffer	pH: 6.50 Preservative: 0.01% Sodium azide Constituents: 0.878% Sodium chloride, 1% BSA, 0.424% Tripotassium orthophosphate
Purity	Affinity purified
Purification notes	Texas Red™ Sulfonil Chloride (Molecular Weight 625 daltons) Absorption Wavelength: 596 nm Emission Wavelength: 620 nm Fluorochrome/Protein Ratio: 4.9 moles Texas Redä per mole of Rabbit IgG
Conjugation notes	This product was prepared from monospecific antiserum by immunoaffinity chromatography using chicken IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities
Clonality	Polyclonal

Applications

The Abpromise guarantee Our **Abpromise guarantee** covers the use of ab6751 in the following tested applications.
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
Immunomicroscopy		Use at an assay dependent concentration.
Flow Cyt		Use at an assay dependent dilution.
IHC-P		Use at an assay dependent concentration.
IHC-Fr		Use at an assay dependent concentration.
ICC/IF		1/1000 - 1/5000. (PMID 18482827).
ELISA		1/10000 - 1/50000.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors