

Product datasheet

Rabbit Anti-Goat IgG H&L (DyLight® 650) preadsorbed ab96918

[1 Image](#)

Overview

Product name	Rabbit Anti-Goat IgG H&L (DyLight® 650) preadsorbed	
Host species	Rabbit	
Target species	Goat	
Specificity	By immunoelectrophoresis and ELISA this antibody reacts specifically with Goat IgG and with light chains common to other Goat immunoglobulins. No antibody was detected against non immunoglobulin serum proteins. Reduced cross-reactivity to chicken, horse, human, mouse, pig, rabbit and rat was detected.	
Tested applications	Suitable for: IHC-P, ICC/IF, Flow Cyt	
Minimal cross-reactivity	Chicken, Horse, Human, Mouse, Pig, Rat	more details
Conjugation	DyLight® 650. Ex: 654nm, Em: 673nm	

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C.
Storage buffer	pH: 6.8 Preservative: 0.09% Sodium azide Constituents: 0.2% BSA, PBS
Purity	Immunogen affinity purified
Purification notes	Antiserum was cross adsorbed using chicken, horse, human, mouse, pig and rat immunosorbents to remove cross reactive antibodies. This antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to DyLight® 650.
Clonality	Polyclonal
Isotype	IgG
General notes	DyLight® 650 replaces DyLight® 649.

Applications


The Abpromise guarantee

Our **Abpromise guarantee** covers the use of ab96918 in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
IHC-P		1/50 - 1/500.
ICC/IF		1/50 - 1/500.
Flow Cyt		1/50 - 1/200.

Images



Emission spectra of DyLight® fluorochromes available in our catalog.

Line colors represent the approximate visible colors of the wavelength maxima.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors