

Product datasheet

Rabbit Anti-Goat IgM H&L (HRP) ab112862

Overview

| | |
|----------------------------|--|
| Product name | Rabbit Anti-Goat IgM H&L (HRP) |
| Host species | Rabbit |
| Target species | Goat |
| Specificity | By immunoelectrophoresis and ELISA this antibody reacts specifically with Goat IgM. This antibody may cross react with IgM from other species. |
| Tested applications | Suitable for: WB, ELISA, ICC/IF, IHC-P |
| Conjugation | HRP |

Properties

| | |
|-----------------------------|---|
| Form | Liquid |
| Storage instructions | Shipped at 4°C. Store at +4°C. |
| Storage buffer | pH: 6.8 Constituents: 0.2% BSA, 0.05% CMIT/MIT based preservative |
| Purity | Immunogen affinity purified |
| Purification notes | Antiserum was solid phase adsorbed to ensure class specificity. The antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to horseradish peroxidase (HRP). |
| Conjugation notes | Molar enzyme/antibody protein ratio is 4:1. |
| Clonality | Polyclonal |
| Isotype | IgG |

Applications

The Abpromise guarantee Our **Abpromise guarantee** covers the use of ab112862 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| Application | Abreviews | Notes |
|-------------|-----------|---|
| WB | | 1/1000 - 1/30000. Colorimetric: 1/1000 - 1/10000; Chemeiluminescent: 1/1000 - 1/30000. |

| Application | Abreviews | Notes |
|-------------|-----------|--------------------------------|
| ELISA | | 1/10000 - 1/100000. (Primary). |
| ICC/IF | | 1/200 - 1/500. |
| IHC-P | | 1/200 - 1/500. |

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors