

## Product datasheet

# Rabbit Anti-Guinea pig IgG H&L (Alkaline Phosphatase) ab6772

## 1 References

### Overview

Product name	Rabbit Anti-Guinea pig IgG H&L (Alkaline Phosphatase)
Host species	Rabbit
Target species	Guinea pig
Tested applications	<b>Suitable for:</b> Dot blot, ELISA, IHC-P, IHC-Fr, Immunomicroscopy, WB, ICC/IF
Immunogen	guinea pig IgG whole molecule
Conjugation	Alkaline Phosphatase

### Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
Storage buffer	Preservative: 0.1% Sodium azide Constituents: 0.878% Sodium chloride, 1% BSA, 0.788% Tris HCl, 0.0095% Magnesium chloride, 0.0014% Zinc chloride, 50% Glycerol (glycerin, glycerine)
Purity	Affinity purified
Purification notes	Alkaline Phosphatase (Calf Intestine) (Molecular Weight 140,000 daltons).
Conjugation notes	This product was prepared from monospecific antiserum by immunoaffinity chromatography using guinea pig IgG coupled to agarose beads
Clonality	Polyclonal

### Applications

**The Abpromise guarantee** Our **Abpromise guarantee** covers the use of ab6772 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes

Application	Abreviews	Notes
Dot blot		
ELISA		
IHC-P		
IHC-Fr		
Immunomicroscopy		
WB		
ICC/IF		

#### Application notes

suitable for immunoblotting (western or dot blot)  
 ELISA and immunohistochemistry as well as other phosphatase-antibody based enzymatic assays requiring lot-to-lot consistency

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

#### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

#### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors