

# Rabbit Anti-Horse IgG H&L (FITC) ab6917

## 2 References

### Overview

Product name	Rabbit Anti-Horse IgG H&L (FITC)
Host species	Rabbit
Target species	Horse
Tested applications	<b>Suitable for:</b> ICC/IF, IHC-P, Immunomicroscopy, Flow Cyt, IHC-Fr, ELISA
Immunogen	Full length native horse whole (purified).
Conjugation	FITC. Ex: 493nm, Em: 528nm

### Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C.
Storage buffer	pH: 6.50 Preservative: 0.01% Sodium azide Constituents: 0.42% Tripotassium orthophosphate, 0.87% Sodium chloride, 1% BSA
Purity	Affinity purified
Purification notes	This product was prepared from monospecific antiserum by immunoaffinity chromatography using horse IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities.
Conjugation notes	Diaminotriazinylaminofluorescein (DTAF) (Molecular Weight 530 daltons) Absorption Wavelength: 495 nm Emission Wavelength: 528 nm Fluorochrome/Protein Ratio: 2.3 moles DTAF per mole of Rabbit IgG
Clonality	Polyclonal

### Applications

**The Abpromise guarantee** Our **Abpromise guarantee** covers the use of ab6917 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
ICC/IF		
IHC-P		
Immunomicroscopy		
Flow Cyt		
IHC-Fr		
ELISA		

#### Application notes

Flow Cyt: Use at an assay dependant concentration.

IHC-P: Use at an assay dependant concentration.

IHC-Fr: Use at an assay dependant concentration.

IM: Use at an assay dependant concentration.

Not tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

#### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

#### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors