

Rabbit Anti-Human IgM H&L (FITC) ab8506

1 References

Overview

Product name	Rabbit Anti-Human IgM H&L (FITC)
Host species	Rabbit
Target species	Human
Tested applications	Suitable for: ICC/IF
Conjugation	FITC. Ex: 493nm, Em: 528nm

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C.
Storage buffer	pH: 7.50 Constituents: 0.164% Sodium phosphate, 50% Glycerol (glycerin, glycerine), 0.58% Sodium chloride
Purity	Immunogen affinity purified
Purification notes	From the crude polyclonal the crossreactive antibodies were extracted by incubation with Sepharose-bound human IgA and IgG. Specific antibodies were absorbed by incubation with Sepharose-bound human IgM. Specific antibodies were eluted by acidic buffer at pH 2.5 followed by neutralisation and dialysis. After repeated binding with immobilized human IgM a minimum of 65% protein bound.
Conjugation notes	The purified polyclonal was conjugated to fluorescein isothiocyanate isomer I followed by gel-filtration and ion-exchange chromatography to clear unbound conjugate.
Clonality	Polyclonal
General notes	F/P ratio is 4-6 mol Fluorescein per 1 mol of goat IgM. There is low background fluorescence in normal tissue in contrast to cases of autoimmune disease where IgM is deposited.

Applications

The Abpromise guarantee Our **Abpromise guarantee** covers the use of ab8506 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
ICC/IF		

Application notes For immunohistochemical staining of human tissue sections a 1:40 dilution is recommended.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors