

Product datasheet

Rabbit Anti-Mouse IgG H&L (Biotin) preadsorbed ab98784

★☆☆☆☆ **1 Abreviews** [1 Image](#)

Overview

Product name	Rabbit Anti-Mouse IgG H&L (Biotin) preadsorbed
Host species	Rabbit
Target species	Mouse
Tested applications	Suitable for: IHC-P, ICC/IF, ELISA, WB
Minimal cross-reactivity	Human, Rat more details
Conjugation	Biotin

Properties


Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
Storage buffer	pH: 6.8 Preservative: 0.09% Sodium azide Constituents: PBS, 0.2% BSA
Purity	Immunogen affinity purified
Purification notes	This antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to Biotin.
Clonality	Polyclonal
Isotype	IgG

Applications

The Abpromise guarantee Our **Abpromise guarantee** covers the use of ab98784 in the following tested applications.
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
IHC-P		1/500 - 1/5000.
ICC/IF		1/200 - 1/1000.
ELISA		1/20000 - 1/400000. as primary
WB		1/20000 - 1/400000. chemiluminescent detection

Images



Immunohistochemical analysis of dog skin tissue, detecting staining of NQO1 with ab98784.

Tissue was fixed with formaldehyde and blocked with none for ; antigen retrieval was by heat mediation in a citrate buffer (pH 6). Samples were incubated with primary antibody (**ab28947 Mouse monoclonal [A180] to NQO1** at 1/175 in BSA in TBS) for 45 minutes. ab98784 (1/500) was used as the secondary antibody.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Rabbit Anti-Mouse IgG H&L (Biotin) preadsorbed (ab98784)

This image is courtesy of an anonymous Abreview

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors