

Product datasheet

Rabbit Anti-Mouse IgG+IgM+IgA H&L ab8516

1 References

Overview

Product name	Rabbit Anti-Mouse IgG+IgM+IgA H&L
Host species	Rabbit
Target species	Mouse
Specificity	Antibodies were raised in rabbit to mouse IgG, IgM and IgA
Tested applications	Suitable for: ELISA, Dot blot

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
Storage buffer	pH: 7.50 Constituents: 0.164% Sodium phosphate, 50% Glycerol (glycerin, glycerine), 0.58% Sodium chloride
Purity	Immunogen affinity purified
Purification notes	Specific antibodies were absorbed by incubation with Sepharose-bound mouse immunoglobulins. Specific antibodies were eluted by acidic buffer at pH 2.5 followed by neutralisation and dialysis. After repeated binding with immobilized mouse immunoglobulins a minimum of 65% protein bound.
Clonality	Polyclonal

Applications

The Abpromise guarantee Our **Abpromise guarantee** covers the use of ab8516 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
ELISA		1/1000.
Dot blot		1/100.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors