

# Rabbit Anti-Mouse IgM mu chain ab9175

## 1 References

### Overview

<b>Product name</b>	Rabbit Anti-Mouse IgM mu chain
<b>Host species</b>	Rabbit
<b>Target species</b>	Mouse
<b>Specificity</b>	The affinity purified anti IgM mu antibodies only react with IgM by immunodiffusion and Immunoelectrophoretic techniques. A 1 mg/ml solution of the antibody is reacted against normal mouse serum and IgG1, IgG2a, IgG2b, IgG3, IgM and IgA. By ELISA, there is minimal reactivity to IgG.
<b>Tested applications</b>	<b>Suitable for:</b> ELISA, Immunoelectrophoresis, WB
<b>Immunogen</b>	Rabbits were immunized with highly purified normal mouse IgM in Freund's adjuvant.

### Properties

<b>Form</b>	Liquid
<b>Storage instructions</b>	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
<b>Storage buffer</b>	pH: 7.20
<b>Purity</b>	Immunogen affinity purified
<b>Purification notes</b>	The antiserum was solid phase absorbed to make it IgM mu chain specific and the anti mu chain antibodies were isolated by solid phase immunoaffinity chromatography.
<b>Clonality</b>	Polyclonal

### Applications

**The Abpromise guarantee** Our **Abpromise guarantee** covers the use of ab9175 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
ELISA		Use at an assay dependent concentration.
Immunoelectrophoresis		Use at an assay dependent concentration.

Application	Abreviews	Notes
WB		Use at an assay dependent concentration.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

#### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

#### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors