

Product datasheet

Rat monoclonal [23G3] Anti-Mouse IgE epsilon chain (Biotin) ab99585

Overview

Product name	Rat monoclonal [23G3] Anti-Mouse IgE epsilon chain (Biotin)
Host species	Rat
Target species	Mouse
Specificity	Reacts with the epsilon heavy chain of BALB/c mouse IgE as demonstrated by ELISA.
Tested applications	Suitable for: Flow Cyt, ELISA, ICC/IF
Conjugation	Biotin

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle. Store In the Dark.
Storage buffer	pH: 7.4 Preservative: 0.1% Sodium azide Constituent: PBS
Purity	Affinity purified
Clonality	Monoclonal
Clone number	23G3
Isotype	IgG1
Light chain type	kappa

Applications

The Abpromise guarantee Our **Abpromise guarantee** covers the use of ab99585 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes

Application	Abreviews	Notes
Flow Cyt		Use 1µg for 10 ⁶ cells. ab18403 - Rat monoclonal IgG1, is suitable for use as an isotype control with this antibody.
ELISA		1/5000 - 1/10000.
ICC/IF		Use at an assay dependent concentration.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors