

Chromatin: Histones, Nucleosomes, Chromosomes and Genomes
February 9–11 2008, Singapore
www.abcam.com/chromatinsingapore

Thank you for registering to attend the forthcoming Chromatin meeting. Below is some information you will find useful to read through prior to your attendance.

Location

The conference will take place at the:

Pan Pacific Hotel
7 Raffles Boulevard
Marina Square
Singapore 039595
Tel: 65 6336 8111
www.panpacific.com/singapore

Registration

The main registration will take place:

Monday February 9 from 09.00 – 12.00

Scientific sessions will commence shortly after lunch at 13:30.

Upon registration you will be provided with the full program and book of abstracts for the meeting. A provisional program is currently available on the website.

Accommodation

The Pan Pacific has been voted Asia's leading Business Hotel for the past 3 years so we're sure you'll enjoy your stay.

For those who have a bedroom reservation, the core nights of 8, 9, and 10 February are booked. You will receive confirmation of your accommodation (including room type and duration of stay) via separate email.

Check-In

Upon arrival at the hotel you should check-in at the main reception. A credit card may be required to cover incidental extras such as telephone charges. Please note you will not receive an individual reservation confirmation from the hotel as your booking has been managed via Abcam. A member of Abcam staff will be available at the hotel in the very unlikely event that you encounter difficulties.

The hotel will make all reasonable efforts to ensure your room is ready for check-in at your arrival time but please note the official check-in time is 3pm

In room facilities

All rooms include an e-fridge, in-room safe, iron and board, flatscreen LCD television with over 30 cable channels, electronic drapes and air-conditioning. The televisions are complemented with Audio Visual portals, allowing technological connectivity with electronic gadgets.

In line with Singapore's no-smoking legislation, smoking is only permitted on the balcony of certain guest rooms where an ashtray is provided for guests' convenience.

Transfers

The quickest and easiest way to the hotel from the airport is by taxi. It is a 20 minute transfer, depending on traffic and will cost approximately \$15.

Posters

Those delegates presenting posters should have received formal confirmation of acceptance. If you have not done so, please contact events@abcam.com.

Each delegate will be provided with a board measuring 2 metres (height) by 1 metre (width) (suitable for standard A0 size posters) and will be notified upon registration where to display. Adhesive will be provided.

Dress Code

The dress code for the conference is smart casual throughout. We kindly request delegates ensure they wear appropriate attire during scientific sessions.

Internet Access

All bedrooms have internet access.

Money and Local Information

The local currency is the Singapore Dollar, although US Dollars, Australian Dollars, Yen and Pounds Sterling are also accepted at most major shopping centers in Singapore. ATMs are widely available and American Express, Diners Club, MasterCard and Visa are accepted.

Time Zone: GMT + 8

Languages: most Singaporeans are bilingual and speak English

Electricity: 220/240 volts AC, 50Hz

Plug fittings of the three-pin square type are in use. Many hotels have 110-volt outlets

Remember your plug adaptor with you for any electrical appliances!!

Visas

Please ensure you have an up to date passport before booking travel and are aware of any visa requirements which may be in operation for nationals of your country. Please see the World Travel Guide for further information.

Travel & Health Insurance

As before any travel, delegates are strongly advised to ensure they have adequate travel and medical insurance. Delegates will be provided with emergency contact numbers upon registration at the conference.

We very much look forward to welcoming you to Singapore.

Abcam Events Team