

Product datasheet

Anti-11-Ketotestosterone antibody ab54042

Overview

Product name	Anti-11-Ketotestosterone antibody
Description	Rabbit polyclonal to 11-Ketotestosterone
Host species	Rabbit
Tested applications	Suitable for: ELISA
Species reactivity	Reacts with: Mouse, Rat, Rabbit, Goat, Chicken, Human, Pig
Immunogen	11-Ketotestosterone conjugated to BSA

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Storage buffer	Preservative: 0.01% Thimerosal (merthiolate) Constituents: 50% Glycerol, PBS, pH 7.5
Purity	Protein G purified
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab54042** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
ELISA		Use at an assay dependent concentration.

Target

Relevance	While testosterone is the primary androgenic male steroid found in mammals, 11-keto-testosterone is a second key androgenic steroid found in fish, where it has been shown to stimulate spermatogenesis.
------------------	--

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors